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zzhhaaooou:

Immunostaining of Cells Adherent to Coverslips

Procedure:

- 1) Immerse 18 mm² glass coverslips in EtOH. In the tissue culture hood, individually pull out and carefully flame to sterilize. Allow to cool then place in 35 mm dishes, or in 6 well plates. Coat overnight (4°C) with 3 ml/dish or well of 40 µg/ml BMS made up in ddH₂O.
- 2) Remove coating solution and block with 1% BSA for 4 hr at 4°C.

3) Isolate cells and plate at 6×10^5 cells/ml (3 ml/well). Allow cells to adhere for 1 hr in incubator. Gently pull off media containing non-adherent cells and wash once with media. Pull off media and add 4% formaldehyde (25 ml of 16% ampule stock made up to 100 ml in PBS) to dishes or wells. Fix for 30 min at room temperature.

4) Remove coverslips. Those for immunostaining are placed vertically in two ceramic coverslip holders. The remainder can be stored at -70°C .

5) Coverslips in ceramic holders are washed two times in PBS. For washes and incubations, use a 250 ml beaker and a 100 ml volume of wash or reagent. Permeabilize by immersion in PBS-Tween ('PBS-T'; PBS containing 0.1% Tween 20) for 15 min at room temperature. Wash four times over 5 min with PBS.

6) Block by immersion in PBS containing 1% BSA for 60 min at room temperature. Immerse one group of coverslips in preimmune sera and the other in immune sera diluted in PBS/1% BSA. Cover and incubate overnight at 4°C .

7) The next day, wash five times over 40 min in PBS. At this time, can immerse in Pierce 'peroxidase suppressor' for 30 min at room temperature, then wash several times in PBS. Immerse in secondary peroxidase-labeled antibody diluted 1/1,000 in PBS containing 1% normal goat serum. Incubate for 60 min at room temp. Wash five times over 40 min in PBS.

5) Place coverslips flat on Parafilm and add Pierce 'metal enhanced DAB' diluted 1/10 in peroxide buffer. Allow reaction to go for 5 - 15 min, then replace in ceramic holders to wash two times in water. Dehydrate in 80%, 90% and 2x 100% EtOH, then immerse in xylene (in chemical hood; all 1 min each). Place cover slip cell layer down on a glass slide containing a drop of mounting medium and examine in the light microscope.

Reagent:

1. PBS-T (1000 ml)

NaCl 8 gm

KH_2PO_4 0.2 gm

Na_2HPO_4 1.15 gm

KCl 0.2 gm

Tween 20 1 ml

Thimerosal 0.1 gm

zzhhaaooou:

Immunostaining of individual cytokine-producing cells smeared on slides

Suitable for use on single-cell suspensions from peripheral blood, lymphoid tissue or cultured cell-lines.

Sample Preparation and Fixation

1. Harvest cells and wash twice in V-bottomed tubes with cold Wash Buffer by centrifugation ($400 \times g$ for 5 minutes) to remove extracellular proteins, including cytokines.
2. Resuspend to $1-5 \times 10^6$ cells/mL in Wash Buffer.
3. Transfer 10-15 of the cell suspension to each reaction field on the adhesion slide.
4. Allow the cells to adhere electrostatically in a monolayer for 10 minutes at room temperature in the humidified chamber to prevent the cells from drying out.
5. Add 50 of ice-cold Fixation Buffer to each field to fix the cells.

6. Incubate for 20 minutes at 4.
7. Wash three times with Wash Buffer to remove free formaldehyde.
8. Add 25 of 2% fetal bovine serum in Wash Buffer to block unbound surface area on the slide.
9. Incubate for 10 minutes at 37.
 - a. Wash slides three times in Wash Buffer-Saponin. The slides are now ready for staining.
 - b. Alternatively, wash the slides in Wash Buffer and allow the slides to dry. Dried slides can be stored at -20 several months before being stained. Prior to staining the slides should be washed in Wash Buffer-Saponin.

Antibody Incubation

All antibody incubations and washes are performed in Wash Buffer-Saponin to keep the cells permeable for antibodies to penetrate the cell membranes.

Detection using Biotin-labeled antibodies

1. Incubate in Endogenous Peroxidase Blocking Buffer for 30 minutes at room temperature in the dark to block endogenous peroxidase activity in the cells (this step can be omitted if cells are to be stained by fluorochromes or non-peroxidase based enzymatic methods).
2. Block endogenous biotin activity with the Avidin/Biotin blocking kit in a two step procedure for 30 minutes in the presence of saponin, described in steps 3-5.
3. Incubate each cell spot on slides with Avidin for 15 minutes supplemented with saponin (0.1% w/v).
4. Wash each cell spot on slides twice in Wash Buffer-Saponin.
5. Incubate in Biotin for 15 minutes supplemented with saponin (0.1% w/v).
6. Wash each cell spot on slides twice in Wash Buffer-Saponin.
7. Incubate each cell spot on slides for 30 minutes at room temperature with 15 unlabeled or biotinylated cytokine-specific antibodies (0.5-5 /mL) diluted in Wash Buffer-Saponin.
8. Wash slides three times in Wash Buffer-Saponin. Note: If using R&D Systems biotinylated antibody skip steps 9 and 10 and continue.
9. Incubate each cell spot on slides for 30 minutes at room temperature with 15 of a biotinylated secondary antibody (either biotin-donkey anti-goat IgG Fab2 diluted 1:700; or biotin-goat anti mouse IgG1 or IgG2A or IgG2B diluted 1:500) in Wash Buffer-Saponin.
10. Wash slides three times in Wash Buffer-Saponin.

Cytokine-specific staining based on either biotinylated primary antibodies or unlabeled primary antibodies along with biotinylated secondary antibodies can then be developed by techniques based on immunofluorescence or immunoenzymatic methods.

Immunoenzymatic Technique

1. Prepare Vectastain Elite ABC-peroxidase reagent according to manufacturers' instructions and supplement with saponin (0.1% w/v), 30 minutes prior to step 2.
2. Incubate for 30 minutes with Vectastain Elite ABC-peroxidase at room temperature in the dark.
3. Wash slides three times with distilled water prior to substrate incubation.
4. Prepare the substrate according to the the protocol given by the manufacturer without any supplementation of saponin.
5. Incubate for approximately 8 minutes (monitor in the microscope) at room temperature in the

dark. A brown color reaction with distinct morphology is developed with DAB in the peroxidase system.

6. Stop the development of the color reaction by repeated washes in distilled water.

7. Counterstain, if desired, with hematoxylin for 1-5 seconds. Slides are then air dried and coverslipped with Immunohistochemistry Mounting Medium. Slides that need to be kept for long periods in optimal shape should be dehydrated in ethanol and mounted in special mounting medium.

zzhhaaouuu:

Frequently Asked Questions

Question 1

Question: Are background signals caused by the developing substrate?Test: Incubate with substrate alone.Cause of Background: Endogenous enzyme (peroxidase/alkaline phosphatase) activity in sampleRemedy: Block with Endogenous Peroxidase Blocking Buffer (with Levamisol in the alkaline phosphate system).

Question 2

Question: Are background signals caused by endogenous biotin activity in the sample?Test: (Block if necessary as in Question 1) Incubate with Vectastain ABC (or ExtraAvidin HRP/AP) and incubate with substrate.Cause of background: Endogenous biotin activity.Remedy: Block with Endogenous Biotin/Avidin Blocking-Buffer.

Question 3

Question: Does the secondary biotinylated antibody cause background signals?Test: (Block if necessary as in Question 1 and/or Question 2). Incubate with the following three steps: secondary antibody; Vectastain ABC-kit (or ExtrAvidin HRP/AP); and substrate.Cause of background: Fc-receptor binding of secondary antibody used at too high concentration.Remedy: Block with species matched serum to secondary antibody or try to dilute the secondary antibody.

Question 4

Question: Does the primary cytokine specific antibody generate background signals?Test: (Block if necessary as in Question 1 and/or Question 2). Use secondary antibody at optimal concentration. Incubate with the following: primary antibody; secondary antibody; Vectastain ABC (or ExtrAvidin HRP/AP); and with substrate.Cause of background: Cytokine-detection antibody used at too high concentration or unsuitable for immunostaining.Remedy: Titrate antibody or try another cytokine detecting antibody.

Question 5

Question: Is the cytokine staining specific?Test: Incubate cytokine-detecting antibody with target cytokine overnight and finally add 0.1 (saponin prior to staining as in Question 4.Cause of Background: Cytokine-detecting antibody used at too high concentration or unsuitable for immunostaining.Remedy: Titrate antibody or try another cytokine detecting antibody.

Contents

Technical Hints

Permeabilization: It is crucial that saponin is present during all antibody incubations and washes to make the staining procedure successful. For detection of intracellular cytokines, the

cytokine-specific antibodies must penetrate through the cell surface membrane, the cytosol, the membranes of the endoplasmic reticulum and the Golgi organelle. The detergent, saponin has been shown to intercalate in the membranes to replace cholesterol and to permeabilize cells in a reversible way, maintaining much of the morphology of the membrane structure of the cell.

Fixation: A solution of phosphate-buffered formaldehyde has been found to preserve cell morphology as well as surface and intracellular antigenicity with minuscule cell aggregation and cell loss. Only fixed cells will stand the effects of detergent treatment.

Controls: Evidence for specificity of the cytokine staining should be based on parallel studies of isotype controls, staining with the secondary antibodies alone and an abolishment of immunoreactivity by preabsorption of the cytokine-specific antibody with the corresponding cytokine protein.

Stimulation of PBMNC for cytokine production: The strongest and most diversified cytokine production is seen after PMA-ionomycin activation of Ficoll separated PBMNC. Both monocytes and lymphocytes will be activated. The cells are co-cultivated with PMA (1 ng/mL), ionomycin (500 nM/mL) and harvested after 4 hours. After 4 hours, one can expect to see IL-1a, IL-1b, IL-1ra, IL-2, IL-3, for various periods. If cells are harvested after 3 hours, one can find monokines such as IL-1a, IL-1b, IL-1ra, IL-6, IL-8, MIP-1a, MIP-1b, and TNF-a. However, cytokines produced by lymphocytes, IL-12 or substantial IL-10 production is not found in the cultures.

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